



1. Blocking
 - 1.1 Add **Blocking Buffer** (1x) to the entire sample area (approximately 40 μ l for each 1cm² area).
 - 1.2 Incubate for 60 min at +37 °C in a pre-heated humidity chamber.

 2. Primary antibody incubation
 - 2.1 Use the provided **Primary Antibody Diluent** (1x) to dilute your primary antibodies.
 - 2.2 Decant the Blocking Buffer and add enough of your antibodies to cover the sample area.
 - 2.3 Incubate for 60 min at +37 °C or overnight at +4 °C in a humidity chamber.
 - 2.4 Decant the antibody solution and wash slides for 3x5 min with 1x TBS-T* in a staining jar under gentle agitation.

 3. Probe incubation
 - 3.1 Prepare the probes by diluting **Probe M1** and **Probe R2** in Probe Diluent (dilute 1:40 each).
 - 3.2 Add enough of the probes to cover the sample area.
 - 3.3 Incubate for 60 min at +37 °C in a pre-heated humidity chamber.
 - 3.4 Decant the solution and wash slides for 3x5 min with 1x TBS-T in a staining jar under gentle agitation.

 4. Reaction A
 - 4.1 Start preparing Reaction A by diluting **Buffer A** (5x) 1:5 in water. Vortex and spin down.
 - 4.2 Add **Enzyme A** (dilute 1:40) Mix gently by pipetting and spin down.
 - 4.3 Add enough Reaction A to cover the sample area.
 - 4.4 Incubate for 60 min at +37 °C in a pre-heated humidity chamber.
 - 4.5 Decant the solution and wash slides for 2x3 min with 1x TBS-T in a staining jar under gentle agitation.

 5. Reaction B
 - 5.1 Start preparing **Reaction B** by diluting Buffer B (5x) 1:5 in water. Vortex and spin down.
 - 5.2 Add **Enzyme B** (dilute 1:40). Mix gently by pipetting and spin down.
 - 5.3 Add enough Reaction B to cover the sample area.
 - 5.4 Incubate for 30 min at 37 °C in a pre-heated humidity chamber.
 - 5.5 Wash slides for 2x3 min with 1x TBS-T in a staining jar under gentle agitation.

 6. Reaction C
 - 6.1 Start preparing **Reaction C** by diluting Buffer C (5x) 1:5 in water. Vortex and spin down.
 - 6.2 Add **Enzyme C** (dilute 1:40). Mix gently by pipetting and spin down.
 - 6.3 Add enough Reaction C to cover the sample area.
 - 6.4 Incubate for 90 min at +37 °C in a pre-heated humidity chamber.
 - 6.5 Decant the solution and wash slides for 2x10 min with 1x TBS in a staining jar under gentle agitation.
 - 6.6 Wash slides for 15 min with 0.1x TBS in a staining jar under gentle agitation.

 7. Mounting
 - 7.1 Tap off excess wash buffer from the slides.
 - 7.2 Mount the slides with a coverslip using water-based anti-fade mounting medium with DAPI.
 - 7.3 Wait for 15 minutes.
 - 7.4 Image your slides in a fluorescence or confocal microscope, using 20x objective or higher and filter set for DAPI and Texas red.
- Protect from light
- Protect from light



KIT COMPONENTS

Box 1: Catalog number ND0110

Material	Product Number	Amount	Storage
Blocking Buffer / Primary Antibody Diluent (1x)	ND00.1.1	8 ml	2–8°C
Probe Diluent (1x)	ND001.2	4 ml	DO NOT FREEZE
Probe M1 (40x)	ND00.1.4	100 µl	
Probe R2 (40x)	ND00.1.3	100 µl	

Box 2: Catalog number ND0120

Material	Product Number	Amount	Storage
Buffer A (5x)	ND00.2.1	800 µl	-15 – -20°C
Enzyme A1 (40x)	ND00.2.2	100 µl	
Buffer B (5x)	ND00.2.3	800 µl	PROTECT FROM LIGHT
Enzyme B (40x)	ND00.2.4	100 µl	
Buffer C (5x)	ND00.2.5	800 µl	
Enzyme C (40x)	ND00.2.6	100 µl	

*TBS-T (Tris-buffered saline supplemented with 0.05% Tween 20)

General guidelines

- Reaction volume depends on the sample size.
- Use best practices when pipetting to minimize reagents consumption.
- Centrifuge vials before pipetting.
- Completely defrost all buffer mixtures (A, B, C) at room temperature, vortex well before use.
- Vortex and spin-down all enzymes (A, B and C) before use.
- Keep enzymes on ice or on a frozen cold block.
- Wait to add enzymes until immediately prior to adding to the sample.
- Buffer C is a light-sensitive reagent. Always keep it protected from light.
- Remove excess washing buffer from samples before adding reagent.
- Do not allow slides/samples to dry.
- Preheat the humidity chamber before each step.
- Incubation times or temperatures other than those specified may give erroneous results.
- As with any product derived from biological sources, proper handling procedures should be used.
- Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
- Unused solution should be disposed according to local regulation.

IMPORTANT: Appropriate precautions should be taken to avoid antibodies cross-contamination. Cross-contamination is the main source of unspecific background due to the high sensitivity of the assay.

Avoid bulk washing methods when multiple antibodies are used.